

a.max AM-1

Capsule Mixer



- Optimal homogeneous mixing results
- Easy loading and unloading
- Noise and vibration free
- High-speed motor

a.max AM-1 Capsule Mixer

Looking to mix capsules of amalgams, glass ionomer, cements, or other capsuled dental materials? a.max AM-1 capsule mixer is the best choice on the market, with a high speed motor to power the unit at 4,800 oscillations per minute, ensuring consistent and homogeneous mixing results. Through patented mechanisms, operation of AM-1 is stable, quiet, and vibration-free. In addition, two mode settings allow users to save their mixing time preferences for easy and quick operation.

Mixing Arms

Redesigned to be compatible with all types of capsules.

Safety Cap

Paired with a built-in microswitch that prevents operation if the cap is open for maximum safety.



High Speed Motor

Mixes at 4,800 opm.

“Set” Button

“Start / Stop” Button

“Time” Button

Power Switch

Specifications:

- Timer Range 1~30 seconds
- Mixing Mode 2
- Speed 4,800 rpm
- Dimensions 210(L) x 200(W) x 150(H) mm
- Weight 3.1 kg
- Power Input AC 100~240V, 47/63 Hz
- Power Output DC 24V \equiv 2.7A



MONITEX INDUSTRIAL CO., LTD

5F, 70, Sec. 1, Guangfu Rd., Sanchong Dist., New Taipei City 24158, TAIWAN
TEL: + 886 2 29993869 FAX: + 886 2 29993647
Email: monitex@monitex.com.tw
www.monitex.com.tw

