





ART³: A NEW WAY TO BULK FILL

Achieve lifelike restorations with the next generation of clinically-proven Adaptive Response Technology (ART) found in Harmonize™ and Optishade Universal Composite. Powered by Enhanced Adaptive Response Technology, SimpliShade Bulk Fill Flow offers easy handling and low shrinkage.



Adaptive Rheological Response Technology for Ease of Sculpting


ART³ uses rheological modifiers that reduce the composite's viscosity as you manipulate it. When you're done shaping, viscosity recovers and anatomy stabilizes before light-curing — reducing the need for resculpting.



Human Enamel SimpliShade Bulk Fill Flow

Adaptive Color-Blending Response Technology for Beautiful Blending

Achieve exceptional blending with special nanofillers that mimic enamel structures and interact with light in a natural way. The restoration and tooth become virtually indistinguishable.



Adaptive Shrinkage Stress Response Technology for Durable Results

Increase confidence in long-lasting restorations. ART³ relieves interfacial stress between the filler and resin matrix during curing. Less stress on the adhesion layer allows for enhanced marginal adaptation.

KEEP IT SIMPLE

Choose SimpliShade Bulk Fill or SimpliShade Universal Composite for your capping layer.

Ordering Information



SimpliShade Bulk Fill Flow

Part #	Item Description
37217	SimpliShade Bulk Fill Flow Syringe
37218	SimpliShade Bulk Fill Flow Unidose 20 Pack



SimpliShade Bulk Fill

Part #	Item Description
37211	SimpliShade Bulk Fill Composite Syringe
37212	SimpliShade Bulk Fill Composite Unidose 20 Pack



SimpliShade Universal Composite

Part #	Item Description
37001	SimpliShade Universal Composite Syringe, Light
37002	SimpliShade Universal Composite Syringe, Medium
37003	SimpliShade Universal Composite Syringe, Dark
37004	SimpliShade Universal Composite Syringe, Bleach White
37005	SimpliShade Universal Composite Syringe, Universal Opaque
37006	SimpliShade Universal Composite Unidose 20 Pack, Light
37007	SimpliShade Universal Composite Unidose 20 Pack, Medium
37008	SimpliShade Universal Composite Unidose 10 Pack, Dark
37009	SimpliShade Universal Composite Unidose 10 Pack, Bleach White
37010	SimpliShade Universal Composite Unidose 10 Pack, Universal Opaque



**One Shade. One Step.
A Better Way to Bulk Fill.**

For more information scan the QR Code or visit
<https://go.kerrdental.com/simplishadebulkfill-jada-0>

SimpliShade™ BULK FILL FLOW



Single-Shade Simplicity Meets Superior Adaptability

Reach every area of the tooth's anatomy while saving time with SimpliShade Bulk Fill Flow. SimpliShade Bulk Fill Flow easily adapts to the cavity floor — with exceptional durability and depth of cure.

kerrdental.com/en-au/
Kerr Australia Pty Limited | Level 4, 7 Eden Park Drive, Macquarie Park | Australia
Australia: 1800 643 603 | New Zealand: 0800 167 626 | orders.anz@kerrdental.com

©2024 Kerr Corporation. SimpliShade and Harmonize are trademarks of Kerr Corporation. All other trademarks are properties of their respective owners.
Data on File. MKT-22-0880-ROW Rev 0



ONE SHADE IS ALL IT TAKES

Get the adaptability, durability, depth of cure, and simplicity you need with SimpliShade Bulk Fill Flow.

Flowable Features Single-Shade Simplicity



One-Shade System

Using one shade that blends seamlessly with all 16 VITA® classical shades, SimpliShade Bulk Fill Flow saves you time selecting shades and managing inventory.



Save Time

Simplify your workflow with just one step to cure up to 4mm, before adding a capping layer.



Go with the Flow

Fill every area of the tooth's anatomy with a self-leveling bulk fill flowable that easily adapts to the cavity floor.



Lower Shrinkage

Helps assure a long-lasting restoration with lower polymerization shrinkage and shrinkage stress when compared to other bulk fill brands.




Add Confidence


Restore confidently with a flowable bulk fill composite that offers comparable or higher mechanical strength than the competitors.

Shrinkage Stress

SimpliShade Bulk Fill Flow offers lower shrinkage stress than Competitor 1, Competitor 2, Competitor 3, and Competitor 4, resulting in lower marginal degradation and postoperative sensitivity for a long-lasting restoration.



Radiopacity (AI%)




Radiopacity

Tests show that SimpliShade Bulk Fill Flow provides higher radiopacity compared to Competitor 1, Competitor 3, and Competitor 4. Higher radiopacity enables clinicians to easily identify restorations radiographically.


Wear

SimpliShade Bulk Fill Flow offers a lower wear rate at interproximal contacts than Competitor 5. Low wear rate results in a long-lasting restoration.



Mechanical Strength

SimpliShade Bulk Fill Flow has comparable or higher mechanical strength than Competitor 1, Competitor 2, and Competitor 5, resulting in better resistance to chipping and fracture over time.



Adaptation

Easily adaptable, SimpliShade Bulk Fill Flow fills every area of the tooth's anatomy.

