



elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

Seamless Connection. Greater Confidence.

Add confidence to your shaping, with intuitive elements Connect and Apex Connect. Designed to work seamlessly together, interconnected system supports patient safety while simplifying your procedures.

elements™ Connect and Apex Connect

An Interconnected System



Dual Light Indicators

With intuitive light indicators in both motor handpiece and the Apex file clip, you'll have greater visibility for more confident shaping.



Patient Safety

The intuitive, easy-to-see smart light and audible indicators, as well as a large, bright viewing display alert you when you're reaching the apex. May be used independently or as a connected system.



Portable Convenience

Easily move between operatories with this interconnected system; lightweight, ergonomic and cordless design.



Versatile Motions

With Adaptive, Reciprocation and Rotary motions in one motor, elements Connect allows offices to maintain their existing techniques and protocols.

Preset

9 Kerr preset file settings, 1 reciprocating setting and 4 customizable settings.

Continuous Use



Two removable batteries are included to ensure one is always charged so you can complete your procedures safely without interruption.




Built-in Accuracy

Ensure measurement accuracy with Apex Connect's automatic self-calibration feature that activates every time the device is turned on.

Shape with Confidence

View real-time working length information during procedures with a smart light indicator that changes color, as well as audible indicators, as you reach the apex.



Ordering Information

elements Connect

815-1901-AU	elements Connect Cordless Motor Kit
815-1906	elements Connect AC/DC Adapter
815-1908	elements Connect Contra Angle
815-1909	elements Connect Charging Base
815-1910	elements Connect Battery Pack

Apex Connect

815-1902-AU	Apex Connect Kit
815-1907	Apex Connect AC/DC Adapter
815-1912	Apex Connect Probe Cord
815-1913	Apex Connect Lip Hook
815-1914	Apex Connect File Holder Adapter
815-1915	Apex Connect File Holder A
815-1916	Apex Connect File Holder B
815-1917	Apex Connect Replacement Battery

Power Cord

815-1904-AU	elements/Apex Connect Power Cord
-------------	----------------------------------



For more information scan the QR Code
or visit to.kerrdental.com/en/elements-connect-apex-connect

**To Learn More About
Kerr Endodontic Solutions
Contact Your Kerr Rep**

