

Kerr™

OptiBond™ Family of Bonding Agents

The Gold Standard for Adhesion



Powered by the GPDM Monomer and Backed with Over 20 Years of Clinical Success.



OptiBond™ Universal

Kerr Part #	Description
36517	Bottle Kit
36519	Bottle Refill
36518	Unidose



OptiBond™ Solo Plus

Kerr Part #	Description
29669	Unidose (pkg of 100 - 10 ml total)
29682	Unidose Kit
31513	Bottle (5 ml)
31514	Bottle Kit



OptiBond™ eXTRa Universal

Kerr Part #	Description
36658	Bottle Kit
36659	Unidose Kit
36660	Bottle Refill - Primer
36661	Bottle Refill - Adhesive



OptiBond™ FL

Kerr Part #	Description
33352	Unidose Kit
26684	Bottle Kit
35265	Bottle Refill - Primer
35266	Bottle Refill - Adhesive



GPDM Monomer

OptiBond GPDM Monomer ensures more effective and reliable adhesion to both the tooth structure and the resin material, creating a better compatibility with both substrates and an overall bond strength.

Kerr's experience in the market (over 130 years), provides us unique insight and expertise that we would love to share with you and your associates! We're committed to making sure that supply to all our partners is continuous and uninterrupted, and our team of dedicated specialists are at your disposal and standing by to assist.

***RISK FREE TRIAL - BUY 2, GET 2 FREE**

OptiBond™ Selection Guide

LEGEND

- ✓ Recommended
- ✓ Secondary Recommendation
- X Not Recommended
- Not Applicable
- Yes
- No



OptiBond™
FL

OptiBond™
Solo Plus

OptiBond™
eXTRa
Universal

OptiBond™
Universal

Indications

	OptiBond™ FL	OptiBond™ Solo Plus	OptiBond™ eXTRa Universal	OptiBond™ Universal
Direct	✓	✓	✓	✓
Indirect	X ¹	✓	✓	✓
Anterior	✓	✓	✓	✓
Posterior	✓	✓	✓	✓

Substrates

	OptiBond™ FL	OptiBond™ Solo Plus	OptiBond™ eXTRa Universal	OptiBond™ Universal
Dentin/Enamel/Uncut Enamel	✓	✓	✓	✓
All Ceramics (Including Lithium Disilicate Ceramic)	X ¹	✓ ²	✓	✓
Metal Based (Including Zirconia Based & Alumina Based)	X ¹	✓ ³	✓	✓ ³
Inlay/Onlays	X ¹	✓	✓	✓
Core Build Up	✓	✓	✓	✓
Cementation of Posts	X ¹	✓	✓	✓
Composite/Metal Repair	✓	✓	✓	✓
Porcelain Repair	✓	✓	✓	✓
Amalgam Sealing	✓	✓	✓	✓

Auxiliary Products Needed

	OptiBond™ FL	OptiBond™ Solo Plus	OptiBond™ eXTRa Universal	OptiBond™ Universal
Gel Etchant	•	•	✓ ⁴	✓ ⁴
Silane	-	•	• ⁵	• ⁵
Metal Primer	-	• ⁵	• ⁵	• ⁵
Total Dark Cure Mode Capability for Cementation	-	•	• ⁶	• ⁶

Bonding Compatibility

	OptiBond™ FL	OptiBond™ Solo Plus	OptiBond™ eXTRa Universal	OptiBond™ Universal
Compatible with all resin cements & core build up	-	• ³	• ⁷	• ³

1. N/A due to High Film Thickness

2. Use with Silane

3. Use only with NX3/Maxcem Elite Cement

4. Optional

5. Not required

6. With NX3/Maxcem Elite/Maxcem Elite Chroma

7. Compatible with any resin cements.

Scan for
Risk Free Trial
T's & C's



Offer ends 31st December 2024