


OptiDisc[™]

Finish as well as you start

Extra-Coarse - CONTOUR



Coarse/Medium - FINISH



Fine - POLISH



Extra-Fine - HIGH GLOSS POLISH











**Buy 2 x OptiDisc Refills
Get 1 x OptiDisc Refill
&
a Small OR Standard
Mandrel FREE**


RISK FREE TRIAL

4 Colours

4 Stages of Abrasiveness

Disc Type	Indicated for	Anterior Convex/Flat area	Posterior Convex/Flat area	Inter-Proximal area
Extra-Coarse	GROSS REDUCTION AND CONTOURING To adjust primary geometric form.			
Coarse/Medium	FINISHING To remove composite excess, shape occlusal anatomy, lingual fissures, and secondary anatomy.			
Fine	POLISHING To erase surface scratches, and reduce surface roughness below Ra = 0.35µm.			
Extra-Fine	HIGH GLOSS POLISHING To reduce surface roughness to high gloss below Ra = 0.2µm.			

The OptiDisc system offers the choice of four grit level discs for contouring, finishing and high gloss polishing. It is suitable for composite, glass ionomers, amalgams, semiprecious and precious metals. Discs are also available in three different diameter sizes depending on your needs: 9.6mm, 12.6mm and 15.9mm.



**Scan for
Risk Free Trial
T's & C's**



Offer ends 31st December 2024