



Oct-Nov 2024 Dealer Promotion

Unit 4, 130-132 Bayfield Road East, Bayswater North,
Victoria 3153, Australia
1800 220 211
sales@matrixdental.com.au
www.matrixdental.com.au

NEW

DURALAST – STERI K & H FILES IN STERILE BLISTER PACKS

Now
Pre-Sterilised!



REDUCE RISK OF NEEDLESTICK INJURIES
Use straight from the perforated packs



SAVE TIME & MONEY
No need to sterilise before use.



NO RISK OF CROSS CONTAMINATION
Single use files will eliminate the risk of cross contamination



GREATER CUTTING EFFICIENCY
Brand new file each and every time



MADE FROM KOREAN SURGICAL STAINLESS STEEL
Quality ISO Coloured coded sizing

Pre sterilised Stainless Steel Files are available in blister packs - both K and Hedstrom files in lengths 21mm, 25mm and 31mm in sizes 6 to 80. Also available in Nickel Titanium.

"ALL PRICING IS INCLUSIVE OF GST"

DURALAST

Steri File

Endo Files Hedstroem

This instrument is directly grinded. High cutting efficiency, use in filing action. Packs of 6.



Ø ISO	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Code: 21mm	80008S	80010S	80015S	80020S	80025S	80030S	80035S	80040S	80045S	80050S	80055S	80060S	80070S	80080S
Code: 25mm	80108S	80110S	80115S	80120S	80125S	80130S	80135S	80140S	80145S	80150S	80155S	80160S	80170S	80180S
Code: 31mm	80708S	80710S	80715S	80720S	80725S	80730S	80735S	80740S	80745S	80750S	80755S	80760S	80770S	80780S

- 80046S 21 mm 45-80 Asstd
- 80116S 25 mm 15-40 Asstd
- 80146S 25 mm 45-80 Asstd
- 80704S 31 mm 45-80 Asstd

Endo Files K

Blade grinded with a square section then twisted in order to obtain cutting power and flexibility. Helix angle higher than on reamer. Packs of 6.



Ø ISO	6	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Code: 21mm	80306S	80308S	80310S	80315S	80320S	80325S	80330S	80335S	80340S	80345S	80350S	80355S	80360S	80370S	80380S
Code: 25mm	80406S	80408S	80410S	80415S	80420S	80425S	80430S	80435S	80440S	80445S	80450S	80455S	80460S	80470S	80480S
Code: 31mm	80606S	80608S	80610S	80615S	80620S	80625S	80630S	80635S	80640S	80645S	80650S	80655S	80660S	80670S	80680S

- 80300S 21 mm 15-40 Asstd
- 80347S 21 mm 45-80 Asstd
- 80417S 25 mm 15-40 Asstd
- 80447S 25 mm 45-80 Asstd
- 80604S 31 mm 45-80 Asstd

\$9.00

Buy 5
Get 1 FREE

DURALAST

ENDODONTIC EXCAVATORS

Each

- Made from spring tempered stainless steel
- Handles 10mm diameter, Satin steel finish, Laser etched markings
- Made in the UK

20%
Discount



27540N

Excavator Endo D/E 31L 1.6mm Duralast Lite

Elongated spoon width 1.6mm. Terminal shank is 15mm long.

\$27.00



27542N

Excavator Endo D/E 33L 2mm Duralast Lite

Elongated spoon width 2.0mm. Curved terminal shank is 15mm long.

\$27.00



27565N

Excavator Endospoon Glick #2 Duralast Lite

Elongated double ended spoon excavator long.

\$30.00

DURALAST
ENDODONTIC PLUGGERS

**20%
Discount**



27916N

Plugger Endo Duralast Lite D/E 5-7

Working ends are 21mm long with a 0.02 taper

- #5 (At 1mm is 0.58mm / At 16mm diameter is 1.2mm)
- #7 (At 1mm diameter is 0.76mm / At 16mm diameter is 1.2mm)

\$18.00



27564N

Plugger Endo Glick #1 Duralast Lite

Plugger has marking at 5mm and 10mm. The paddle is used to place materials.

\$30.00

DURALAST
ENDODONTIC SPREADERS

**20%
Discount**



27560N

Spreader Endo Duralast Lite S/E D11 0.30mm

The working end is 21mm long. At 1mm diameter is 0.30mm. At 16mm the diameter is 0.91mm.

\$18.00



27561N

Spreader Endo Duralast Lite S/E D11T 0.25mm

The working end is 21mm long. At 1mm diameter is 0.25mm. At 16mm the diameter is 0.86mm.

\$20.00

MGuard

LATEX FREE DAM

15cmx15cm Lt/Blue Med

68100

Box-20

- Non latex dental dam made from polyisoprene
- Extremely flexible and doesn't tear
- Easy to place and floss through without the need for lubrication
- Available in Medium Grade Light Blue



**Buy 5
Get 1 FREE**

\$38.50 Box-20

Or **\$35.00** When you Buy +10

ACERO S3

CROWN STAINLESS STEEL

- Stainless Steel Paediatric molar crowns
- Made from Surgical Grade 316L Stainless steel
- Shallow occlusal anatomy
- Universal numbering system
- Supplied in snap lock transportation Kit
- Made in the USA

22000 3S Primary Kit-240

\$1349.00

**40% OFF
Per Kit**

22001 3S Primary Kit-48

\$319.00

**30% OFF
Per Kit**



**Special
Kit Price**

FRASACO

CROWN FORMS STRIP

- Ultra thin temporary celluloid crowns
- Easy to trim
- Available in 58 adult and paedodontic moulds and in various sizes

FR-STKK84 Pedo Strip Crown Form Kit-84

\$275.00

FR-STFBK58 Adult Ant/Post Strip Crown Form Kit-58


\$200.00

FR-STFK174 Adult Anterior Strip Crown Form Kit-174

\$520.00

FR-STBK108 Adult Post Asst Strip Crown Form Kit-108

\$335.00





DURALAST

COMPOSITE TITANIUM NITRIDE COATED

- Titanium Nitride Coated tips for ultra smooth non stick surface
- Handles 10mm diameter, satin steel finish, laser etched markings



REPF14T Composite Tit/Nitride Coat PFI 4T DE#4

Identical flared paddles 11mm long/2.0mm wide

\$60.00



27704N Composite Tit/Nitride Coat PFI #6 Long Paddle/Plug

Double ended paddle 8.4mm x 13.7mm & 6.9mm x 10mm

\$38.00



27708N Composite Tit/Nitride Coat PFI #18 Ball Lite HNDL

Double Ended ball burnisher 1.6mm/3.4mm

\$38.00



REPF139T Composite Tit/Nitride Coat PFI 139T Placement

Identical Flat medium paddles at right angles to each other 9mm long/2.0mm wide

\$38.00

ABSORBENTS

Reflective Shields Plus

- More absorbent than the leading brands
- Does not swell for superior accessibility in small mouths
- Comfortable for the patient
- Unique reflective film allows illumination
- Less expensive than some brands
- Small/Adult size

600700

Reflective Shields Plus

Absorbent Pad

Adult

\$27.00

600710

Reflective Shields Plus

Absorbent Pad

Small

\$27.00



FLOSSERS

REACH ACCESS FLOSSER

- Reusable Flosser: The Reach Access Snap-on Flosser is a reusable flosser
- Refreshing mint flavour

33040

Flosser Cleanpaste Starter C42031036
Box-6



\$45.00

- Reach Starter Flosser comes with refill heads
- Ideal for patients that find flossing difficult
- Refill heads available

33041

Flosser Cleanpaste Refill C42031037
Box-6



\$63.00

**Buy 7
Get 1 FREE**

HOGENSPLITZE

NEEDLES

Box-100


- Tri-Bevel cut for patient comfort
- Siliconised for smooth insertion
- Hub marker shows needle bevel
- Two gauges and 7 lengths to choose from

- Z1015 Needles Anaes 27G S 21mm Imp Thread
- Z1023 Needles Anaes 27G S 25mm Imp Thread
- Z1021 Needles Anaes 27G L 35mm Imp Thread
- Z1019 Needles Anaes 27G L 40mm Imp Thread
- Z1020 Needles Anaes 30G L 25mm Imp Thread
- Z1016 Needles Anaes 30G S 21mm Imp Thread
- Z1022 Needles Anaes 30G XS 12mm Imp Thread



**Can be mixed
for Bulk
Purchase
Price**

\$15.00 Each
or **\$14.00** When You Buy +10
or **\$12.00** When You Buy +20




COREFLON

PTFE Non-Absorbable Suture

Box

- . PTFE Suture swagged to black needles
- . Outstanding visibility, even under high illumination
- . Round body taper cut needles for the smallest of tissue puncture
- . White PTFE monofilament material reduces bacterial tracking -less infections
- . Smooth PTFE is more comfortable for patients
- . Good knot security



- 19000 PTFE 3/0 1/2 Circle Taper Cut 16mm Black Needle
- 19001 PTFE 3/0 1/2 Circle Taper Cut 19mm Black Needle
- 19002 PTFE 4/0 1/2 Circle Taper Cut 16mm Black Needle
- 19003 PTFE 4/0 1/2 Circle Taper Cut 19mm Black Needle
- 19004 PTFE 5/0 1/2 Circle Taper Cut 16mm Black Needle
- 19005 PTFE 6/0 1/2 Circle Taper Cut 13mm Black Needle

Needle	Length	Shape	3/0	4/0	5/0	6/0
DS 1/2 circle Taper Cut	13mm					19005
	16mm		19000	19002	19004	
	19mm		19001	19003		

Buy 1 Box **\$285.00**
Or Buy 3 Box **\$255.00**
Or Buy 10 Box **\$225.00**

20%
Discount

DURALAST

SURGICAL

Each

- Working tips made from 420 grade stainless steel
- High strength and corrosion resistant stainless steel
- Handles are made from Resin and colour coded for ease of identification
- Available in 8 sizes
- Sterilise by Steam only
- Made in the UK

27431R

Elevator Luxation 2c Cvd 2mm Resin Handle



27430R

Elevator Luxation 2s Str 2mm Resin Handle



27433R

Elevator Luxation 3c Cvd 3mm Resin Handle



27434R

Elevator Luxation 3s Str 3mm Resin Handle



27435R

Elevator Luxation 5c Cvd 5mm Resin Handle



\$46.00EA

**20%
Discount**

DURALAST

SURGICAL

Each

- Working tips made from 420 grade stainless steel
- High strength and corrosion resistant
- Handles are made of stainless steel
- Available in 4 sizes
- Made in the UK


27433

Elevator Luxation 3c Cvd 3mm



27434

Elevator Luxation 3s Straight 3mm



27435

Elevator Luxation 5c Cvd 5mm



\$46.00EA

**25%
Discount**

LEDERMIX

- . Over 40 years history using Triamcinolone acetonide & Demeclocycline
- . Fast relief from pain and effectively reduces microbial contamination
- . Used for acute irreversible pulpitis, deep carious lesions infected pulp & necrosis of pulp

Ledermix Paste 5G
007698



\$140.00EA

**15%
Discount**

LEDERMIX

- . Hard setting material containing Triamcinolone acetonide and Eugenol
- . Used for teeth with reversible pulpitis as sedative dressing or lining under filling or to reduce inflammation under temporary crown
- . Used in primary teeth as root filling material

Ledermix Refill 3
Cement Powder 3g
007699

\$94.00EA

Ledermix Hardner
N 2.5ml
007696

\$14.00EA

